

TRANSMISSION OF WRATTEN FILTERS

Compiled by Allie C. Peed, Jr. for The Eastman Kodak Company

Data condensed from Kodak Wratten Filters for Scientific and Technical Use published by the Eastman Kodak Company, manufacturers of the filters.

The following pages give (1) percentage luminous transmittance at wave lengths from 400 to 700 μ at intervals of 10 μ for the standard illuminant "C" adopted by the International Commission of Illumination, (2) dominant wavelength in millimicrons, and (3) percentage of excitation purity. Values of wave length followed by "c" indicate the complementary wave lengths of purple filters which do not have a dominant wave length.

All colorimetric specifications are based on the 1931 standard ICI colorimetric and luminosity data.

The transmittance data are given as representing standard samples of the filters. They are intended only for the information of users in choosing filters which will meet their requirements. Values taken from the tables of data should not be used by research workers as representing precisely the absorption characteristics of a particular filter. If such precise data are needed, they should be determined for the particular filter being used.

Where the spectra extend into the ultraviolet this fact is indicated by an asterisk (*) in the transmission tables immediately beneath the filter number, and quantitative data are not given. The manufacturer should be consulted for this information. Transmission in the ultraviolet of wave lengths less than 330 μ will be eliminated in the case of cemented filters, as glass absorbs ultraviolet radiation of wave lengths shorter than about 330 μ .

Stability ratings are given as three letter combinations following the filter description in the table below. In establishing the stability classifications each filter is exposed to a selected light source for a specific time interval. The following grading system is used to describe the result:

- Class A — stable
- Class B — relatively stable
- Class C — somewhat unstable
- Class D — unstable

The classification letters, for example, AAA, describe the stability to the following three exposure tests in this order:

1. Two weeks' exposure to daylight in a south window
2. Twenty-four hours' exposure to a "Fade-Ometer"
3. Two weeks' exposure at two feet from a 1000-watt tungsten lamp.

Filters are supplied in two forms: as lacquered gelatin film, or as a gelatin film cemented between pieces of optical glass. Filters in glass are cemented between sheets of plane-parallel glass, which is surfaced in quantities and is of sufficient accuracy for general photographic work, and for most scientific purposes.

Most Wratten Gelatin Filters are stocked in 2- or 3-inch squares. Stocks of 2- or 3-inch square filters cemented in glass are maintained only in filters usually used for general photographic work.

The booklet "Kodak Filters and Lens Attachments" gives more valuable information on this subject.

FILTER DATA

No.	Description, use, and stability	No.	Description, use, and stability
	Colorless	25	A — Tricolor red for direct color separation. Contrast effects in commercial photography and in outdoor scenes. Two-color general viewing. Aerial infrared photography and haze cutting, AAA.
0	For compensating thickness of other gelatin filters in optical systems, AAA.	26	Stereo red, AAA.
1	Absorbs ultraviolet below 360 m μ , DDD.	29	Red color separation from transparencies and for the Kodak Fluorescence Process. Strong contrast effects. Copying blueprints. Tungsten tricolor projection, AAA.
1A	Kodak Skylight Filter — Reduces excess bluishness in outdoor color photographs in open shade under a clear, blue sky, ACA.		Magentas and Violets
	Yellow	30	Green absorption, BBC.
2B	Absorbs ultraviolet below 410 m μ , ACA.	31	Green absorption, CCA.
3	Light yellow, CCD.	32	Minus green, CCD.
3N5	No. 3 plus 0.5 neutral density, AAA.	33	Strong green absorption, CCB.
4	Light yellow — Approximate correction on panchromatic materials for outdoor scenes, including sky, CCC.	34	Violet, CDD.
6	K1 — Light yellow — Partial correction outdoors, BBA.	34A	Blue separation — Kodak Fluorescence Process, DCC.
8	K2 — Yellow — Full correction outdoors on Type B panchromatic materials. Widely used for proper sky, cloud, and foliage rendering. Green separation for Fluorescence Process, AAA.	35	Contrast in microscopy, CDD.
	Oranges and Reds	36	Dark violet, CCC.
8N5	No. 8 plus 0.5 neutral density, AAA.		Blues and Blue-greens
9	K3 — Deep yellow. Moderate contrast in outdoor photography (with black-and-white films), AAA.	38	Red absorption, BCA.
11	X1 — Greenish yellow. Correction for tungsten light on Type B panchromatic materials; also for daylight correction with Type C panchromatic materials in making outdoor portraits, darkening skies, or lightening foliage, AAA.	38A	Red absorption. Increasing contrast in visual microscopy, BBB.
12	Minus blue. Haze cutting in aerial photography, AAA.	39	Contrast control in printing motion-picture duplicates (glass) AAA.
13	X2 — Yellow green. Correction for Type C panchromatic materials in tungsten light, ABA.	40	Green for two-color photography (tungsten), CBC.
15	G — Deep yellow. Overcorrection in landscape photography. Contrast control in copying and in aerial infrared photography, AAA.	44	Minus red — Two-color general viewing, DDD.
16	Blue absorption, AAB.	44A	Minus red, DDD.
18A	Transmits ultraviolet and infrared only (glass), AAA.	45	Contrast in microscopy, DDD.
	Oranges and Reds	45A	Highest resolving power in visual microscopy, CDC.
21	Blue and blue-green absorption, CBB.	46	Blue projection (experimental), DDD.
22	Yellow-orange. For increasing contrast in blue preparations in microscopy. Mercury yellow, BAC.	47	Tricolor blue for direct color separation and from Kodak Ektacolor Film for Dye Transfer. Contrast effects in commercial photography. Tungsten and white-flame-arc tricolor projection, BBC.
23A	Light red. Two-color projection — contrast effects, BAB.	47B	Tricolor blue for color separation from transparencies and from Kodak Ektacolor Film for Graphic Arts, BBB.
24	Red for two-color photography (daylight or tungsten). White-flame-arc tricolor projection, AAB.	48	Green and red absorption, CBC.
	Greens	48A	Green and red absorption, AAB.
		49	Dark blue, CCB.
		49B	Very dark blue, BBB.
		50	Very dark blue. Mercury violet, CCC.
		52	Light green, AAB.

TRANSMISSION OF WRATTEN FILTERS (Continued)

FILTER DATA

No.	Description, use, and stability	No.	Description, use, and stability
53	Medium green, CCB Very dark green, AAA.	†86B	Yellowish. Photometric filter (visual), BCA.
55	Stereo green, BBC.	†86C	Yellowish. Photometric filter (visual), AAA.
56	Very light green, CBC.		Light Balancing
57	Green for two-color photography (daylight), CBC.	80A	For Kodachrome Film, Daylight Type, and photographic flood lamps, ABA.
57A	Light green, BBC.	81	Yellowish. For warmer color rendering.
58	Tricolor green for direct color separation. Contrast effects in commercial photography and microscopy, BBC.	81A	Yellowish. For Kodak Ektachrome Film, Type B, with photographic flood lamps.
59	Green for tricolor projection (white-flame-arc), BBB.	81B	Yellowish. For warmer color rendering.
59A	Very light green, BBB.	81C	Yellowish. For Kodachrome Film, Type A, with flash lamps.
60	Green for two-color photography (tungsten), BDC.	81D	Yellowish. For Kodachrome Film, Type A, with flash lamps.
61	Green color separation from transparencies and Kodak Ektacolor Film. Tricolor projection (tungsten), ABC.	81EF	Yellowish. For Kodak Ektachrome Film, Type B, with flash lamps.
64	Red absorption (light), CDB.	82	Bluish. For cooler color rendering.
65	Red absorption, ADB.	82A	Bluish. For Kodachrome Film, Type A, with 3200 K lamps.
65A	Red absorption, CCD.	82B	Bluish. For cooler color rendering.
66	Contrast effects in microscopy and medical photography, DDC.	82C	Bluish. For cooler color rendering.
67A	Red absorption (light). Two-color projection, CDC. Narrow-band	83	Yellowish. For 16 mm Commercial Kodachrome Film and daylight exposure, BBB.
70	Dark red. Infrared photography. Color separation for Kodak Ektacolor Film (with tungsten), ABC.	85	Orange. For Type A Kodak color films and daylight exposure, BAA.
72B	Dark orange-yellow, CCC.	85B	Orange. For Kodak Ektachrome Film, Type B, and daylight exposure, BAB.
73	Dark yellow-green, ABB.		Miscellaneous
74	Dark green. Mercury green, BBC.	79	Photographic sensitometry. Corrects 2360 K to 5500 K, AAA.
75	Dark blue-green, ACC.	87	For infrared photography. Absorbs visual.
76	Dark violet (compound filter), DDD.	87C	Absorbs visual, transmits infrared.
77	Transmits 546 m μ mercury line (glass plus gelatin), AAA.	88A	For infrared photography. Absorbs visual.
77A	Transmits 546 m μ mercury line (glass plus gelatin), AAA.	89B	For infrared photography, AAA.
	Photometrics	90	Narrow-band viewing filter for judging brightness scale of scenes, CCD.
78	Bluish. Photometric filter (visual), BAB.	96	Neutral filters for controlling luminance, AAB.
78AA	Bluish. Photometric filter (visual), BAA.	97	Dichroic absorption, AAA.
78A	Bluish. Photometric filter (visual), AAA.	102	Correction filter for Barrier-layer photocell, ABA.
78B	Bluish. Photometric filter (visual), AAA.	106	Correction filter for S-4 type photocell, AAA.
78C	Bluish. Photometric filter (visual), BAA.		
86	Yellowish. Photometric filter (visual), BBA.		
86A	Yellowish. Photometric filter (visual), AAA.		

Percent transmittance

Wave length	No. 0	No. 1	No. 1A	No. 2B	No. 3	No. 3N5	No. 4	No. 6	No. 8	No. 8N5	No. 9	No. 11	No. 12
400	88.0	85.0	59.0	19.0	—	—	—	7.40	—	—	—	—	—
10	88.5	85.5	76.0	48.0	—	—	—	8.32	—	—	—	—	—
20	88.9	86.0	82.0	67.0	—	—	—	10.4	—	—	—	0.16	—
30	89.3	86.5	84.6	75.3	0.36	—	—	13.5	—	—	—	0.29	—
40	89.6	87.0	86.0	80.0	1.78	—	—	18.9	—	—	—	0.56	—
50	89.8	87.4	86.8	83.0	11.5	1.59	—	27.6	—	—	—	1.32	—
60	89.9	87.8	87.2	85.2	38.0	9.40	6.9	39.0	0.25	0.16	—	4.00	—
70	90.1	88.2	87.5	86.7	68.0	18.5	42.0	52.3	5.50	2.0	1.78	12.0	—
80	90.3	88.5	87.3	88.1	80.8	23.5	74.0	65.8	19.0	6.3	8.31	26.0	—
90	90.4	88.7	86.8	88.8	85.2	25.5	84.7	76.8	41.0	13.2	20.7	43.7	—
500	90.5	89.8	86.3	89.5	86.9	26.3	87.5	83.5	63.5	20.0	34.5	55.0	1.50
10	90.6	89.1	85.5	89.9	87.8	26.7	88.5	87.0	78.0	24.3	48.8	60.0	17.3
20	90.7	89.3	84.8	90.3	88.4	27.0	89.1	88.4	84.1	26.7	62.0	60.2	55.0
30	90.7	89.5	84.3	90.5	89.0	27.2	89.4	89.0	86.5	28.0	76.0	57.8	77.8
40	90.8	89.7	84.0	90.6	89.5	27.5	89.6	89.6	87.7	28.6	83.8	54.2	86.0
50	90.8	89.9	83.9	90.7	89.8	27.8	89.8	89.7	88.4	29.0	87.0	50.0	88.4
60	90.9	90.1	84.1	90.8	90.1	27.9	89.9	89.9	88.8	29.3	88.3	44.8	89.4
70	90.9	90.2	84.8	90.9	90.4	28.0	90.2	90.1	89.2	29.5	88.8	38.9	89.7
80	90.9	90.3	86.0	90.9	90.6	28.4	90.4	90.3	89.5	29.6	89.1	33.1	90.1
90	91.0	90.4	87.4	91.0	90.7	29.0	90.6	90.5	89.8	29.8	89.3	27.6	90.3
600	91.0	90.5	88.5	91.1	90.8	29.5	90.8	90.6	90.1	29.9	89.5	22.7	90.4
10	91.0	90.5	89.5	91.2	90.9	29.5	90.9	90.7	90.3	29.6	89.7	19.0	90.5
20	91.0	90.6	90.2	91.3	91.0	29.3	91.0	90.8	90.5	29.4	89.8	14.9	90.7
30	91.0	90.6	90.6	91.3	91.0	29.1	91.1	90.9	90.7	29.1	89.9	11.4	90.8
40	91.1	90.7	90.8	91.4	91.1	29.0	91.2	91.0	90.9	28.8	90.0	9.10	90.9
50	91.1	90.7	91.0	91.4	91.2	20.4	91.3	91.1	91.0	28.9	90.1	8.05	91.0
60	91.1	90.8	91.1	91.5	91.3	29.6	91.4	91.2	91.1	29.2	90.1	7.50	91.1
70	91.1	90.8	91.1	91.5	91.4	29.8	91.5	91.2	91.2	29.4	90.2	7.05	91.2
80	91.1	90.9	91.1	91.6	91.5	30.0	91.5	91.3	91.3	29.5	90.2	6.50	91.2
90	91.1	90.9	91.1	91.7	91.6	30.2	91.6	91.4	91.4	29.7	90.3	6.10	91.2
700	91.1	91.0	91.1	91.8	91.7	31.0	91.6	91.3	91.5	30.2	90.3	6.20	91.3
Luminous transmit.	90.8	89.9	85.9	90.5	88.3	27.4	87.8	87.5	82.7	27.0	76.6	40.2	73.8
Dominant wave lgth.	571.0	575.0	498.0	570.0	569.5	570.5	569.5	570.3	571.8	572.0	574.4	550.3	576.1
Excitation purity.	0.8	1.5	1.2	5.7	50.0	56.3	64.0	44.7	85.2	84.0	91.4	60.7	97.8

* Some transmission below 400 m μ . Consult the manufacturer.

TRANSMISSION OF WRATTEN FILTERS (Continued)

Wave length	Percent transmittance												
	No. 13	No. 15	No. 16	No. 18A	No. 21	No. 22	No. 23A	No. 24	No. 25	No. 26	No. 29	No. 30	No. 31
400	—	—	—	—	—	—	—	—	—	—	—	48.6	13.8
10	—	—	—	—	—	—	—	—	—	—	—	47.4	14.5
20	—	—	—	—	—	—	—	—	—	—	—	48.5	16.4
30	0.18	—	—	—	—	—	—	—	—	—	—	50.1	25.5
40	0.50	—	—	—	—	—	—	—	—	—	—	49.4	42.7
50	1.35	—	—	—	—	—	—	—	—	—	—	43.0	50.2
60	4.08	—	—	—	—	—	—	—	—	—	—	26.5	40.4
70	11.0	—	—	—	—	—	—	—	—	—	—	13.8	22.6
80	23.5	—	—	—	—	—	—	—	—	—	—	5.00	8.20
90	39.0	—	—	—	—	—	—	—	—	—	—	0.63	1.85
500	50.8	—	—	—	—	—	—	—	—	—	—	—	0.12
10	55.2	1.00	—	—	—	—	—	—	—	—	—	—	—
20	56.5	16.0	3.00	—	—	—	—	—	—	—	—	—	—
30	55.0	52.1	22.0	—	—	—	—	—	—	—	—	—	—
40	51.0	70.7	48.0	—	—	—	—	—	—	—	—	—	—
50	46.0	84.3	69.5	—	2.50	—	—	—	—	—	—	—	—
60	39.2	87.5	79.5	—	29.0	0.25	—	—	—	—	—	—	—
70	32.0	88.7	84.0	—	65.0	19.0	—	—	—	—	—	—	—
80	25.1	89.3	86.3	—	80.6	60.0	11.0	—	—	—	—	0.10	—
90	18.2	89.7	87.8	—	85.4	81.0	47.0	4.55	—	—	—	45.0	—
600	13.5	90.0	90.0	—	87.3	87.0	69.6	37.3	12.6	2.90	—	76.0	0.63
10	9.60	90.1	89.0	—	88.1	88.5	82.7	72.3	50.0	30.0	—	87.4	26.0
20	6.40	90.2	90.0	—	88.7	89.0	85.8	82.9	75.0	63.2	10.0	89.5	67.2
30	3.66	90.3	90.2	—	89.0	89.5	87.2	86.4	82.6	78.9	45.3	90.2	84.0
40	2.20	90.4	90.3	—	89.5	89.8	87.9	87.8	85.5	84.0	71.4	90.5	88.1
50	1.58	90.5	90.4	—	89.9	90.0	88.5	88.5	86.7	86.1	82.7	90.7	89.8
60	1.74	90.6	90.5	—	90.4	90.2	89.4	89.0	87.6	87.2	86.6	90.8	90.2
70	2.62	90.6	90.6	—	90.5	90.3	89.6	89.3	88.2	88.1	88.4	90.9	90.4
80	3.55	90.7	90.7	—	90.5	90.4	89.8	89.9	89.0	88.5	89.4	91.0	90.5
90	4.48	90.7	90.8	0.25	90.6	90.5	90.0	90.2	89.3	89.2	90.3	91.1	90.7
700	5.25	90.8	90.8	1.20	90.6	90.6	90.2	90.3	89.5	89.5	90.4	91.1	90.8
Luminous transmit.	34.5	66.2	57.7	0.0014	45.6	35.8	25.0	17.8	14.0	11.7	6.3	26.6	12.9
Dominant wave lgth.	542.0	579.3	582.7	700.0	588.9	595.1	602.7	610.6	615.1	619.0	631.6	498.6	513.1
Excitation purity.	57.5	99.0	99.3	100.0	99.9	99.9	100.0	100.0	100.0	100.0	100.0	62.4	81.9

Wave length	Percent transmittance												
	No. 32	No. 33	No. 34	No. 34A	No. 35	No. 36	No. 38	No. 38A	No. 39	No. 40	No. 44	No. 44A	No. 45
400	38.0	0.85	64.0	—	48.0	36.5	60.5	33.4	85.2	—	0.44	2.52	—
10	37.9	0.71	70.1	0.1	57.0	45.5	66.5	41.2	78.2	—	0.36	3.39	—
20	40.0	1.17	72.0	40.0	57.6	45.5	72.5	53.0	70.5	—	0.63	6.30	—
30	43.0	1.69	68.4	69.7	47.5	32.7	75.3	58.0	63.3	—	3.63	17.4	—
40	55.5	5.36	58.2	68.7	29.5	15.2	76.2	58.8	53.6	—	13.1	32.7	5.00
50	66.0	14.3	42.3	56.2	12.3	3.7	75.9	57.6	42.5	—	25.4	41.8	19.0
60	66.0	12.4	25.2	40.5	3.5	0.35	74.8	55.2	28.5	3.16	36.5	48.1	29.5
70	57.0	5.00	12.1	23.8	0.25	—	73.4	51.9	17.3	21.6	46.5	51.7	34.4
80	40.0	0.50	2.7	9.2	—	—	71.6	48.5	10.2	44.7	53.6	52.9	35.7
90	21.0	—	0.2	2.3	—	—	69.5	44.6	4.00	61.4	56.8	52.2	34.5
500	9.56	—	—	0.33	—	—	66.7	40.2	1.33	70.2	55.8	49.8	29.7
10	2.51	—	—	—	—	—	63.9	35.8	0.35	72.4	50.9	44.8	21.5
20	0.13	—	—	—	—	—	60.8	31.7	—	70.5	42.1	36.8	11.5
30	—	—	—	—	—	—	57.0	27.2	—	64.8	30.5	26.8	3.80
40	—	—	—	—	—	—	52.6	22.3	—	55.5	18.6	16.8	0.85
50	—	—	—	—	—	—	48.0	17.6	—	44.2	8.99	8.20	—
60	—	—	—	—	—	—	42.8	12.9	—	32.5	3.59	2.95	—
70	—	—	—	—	—	—	37.0	8.78	—	20.3	0.80	0.91	—
80	—	—	—	—	—	—	30.6	5.65	—	9.56	—	0.10	—
90	—	—	—	—	—	—	25.5	3.48	—	3.20	—	—	—
600	6.04	—	—	—	—	—	20.9	2.09	—	1.10	—	—	—
10	41.0	0.80	—	0.13	—	—	16.8	1.15	—	0.32	—	—	—
20	75.0	24.9	—	1.0	—	—	12.9	0.59	—	—	—	—	—
30	86.1	60.8	—	6.3	—	—	10.0	0.28	—	—	—	—	—
40	89.0	78.0	0.4	22.0	—	—	7.79	0.13	—	—	—	—	—
50	90.0	85.0	4.0	45.0	0.1	—	6.68	—	—	—	—	—	—
60	90.6	87.5	20.7	65.0	3.0	0.21	6.20	—	—	—	—	—	—
70	90.7	88.7	45.2	77.3	19.0	7.5	5.91	—	—	—	—	—	—
80	90.8	89.4	66.5	85.0	43.5	29.0	5.41	—	0.50	0.80	—	—	—
90	90.9	89.8	78.8	88.2	66.0	55.0	4.90	—	4.06	6.99	0.18	—	—
700	91.0	90.0	85.0	89.8	77.7	71.3	5.00	—	17.8	23.5	1.60	—	1.00
Luminous transmit.	12.5	5.2	1.3	2.9	0.45	0.25	42.5	17.3	1.2	33.6	15.6	14.4	5.2
Dominant wave lgth.	551.7	498.0	424.0	564.8	566.8	566.4	483.5	478.9	450.6	516.2	589.1	483.4	481.5
Excitation purity.	79.6	88.3	94.4	91.4	96.3	97.8	41.8	69.8	98.9	48.5	72.9	77.2	88.4

* Some transmission below 400 m μ . Consult the manufacturer.

TRANSMISSION OF WRATTEN FILTERS (Continued)

Wave length	Percent transmittance												
	No. 45A	No. 46	No. 47	No. 47B	No. 48	No. 48A	No. 49	No. 49B	No. 50	No. 52	No. 53	No. 54	No. 55
400	—	1.20	7.80	16.0	0.96	5.65	3.30	1.70	0.45	2.18	—	—	—
10	—	0.60	17.4	29.5	3.16	10.0	4.28	2.00	0.39	1.51	—	—	—
20	—	0.80	34.0	43.6	8.25	16.0	6.93	3.55	0.59	0.80	—	—	—
30	1.00	5.98	47.0	50.0	15.0	21.0	11.2	7.00	2.63	0.44	—	—	—
40	8.81	19.0	50.3	47.2	22.6	25.0	18.9	13.0	8.90	0.41	—	—	—
50	17.4	30.1	48.3	36.0	30.3	26.2	25.6	17.4	14.0	0.69	—	—	—
60	20.9	33.8	43.4	25.0	33.2	22.9	24.0	14.8	12.3	1.45	—	—	0.20
70	21.6	32.1	36.2	13.2	29.6	16.5	15.7	7.60	5.36	2.70	0.10	—	2.90
80	20.5	27.0	28.5	4.5	22.4	9.55	6.93	2.76	1.55	4.90	0.7	—	13.1
90	18.0	20.2	19.6	1.3	14.1	4.27	2.14	0.40	0.10	8.50	2.14	—	34.2
500	14.4	11.1	0.36	0.17	7.30	1.58	0.46	—	—	13.3	4.47	—	53.4
10	10.1	4.39	—	—	2.64	0.48	—	—	—	18.2	7.24	0.10	67.0
20	5.60	1.66	—	—	0.50	—	—	—	—	23.7	10.7	0.31	69.3
30	2.52	0.35	—	—	—	—	—	—	—	28.5	14.0	0.64	65.1
40	0.4	—	—	—	—	—	—	—	—	32.1	16.6	0.89	56.7
50	0.10	—	—	—	—	—	—	—	—	33.1	17.3	0.93	45.0
60	—	—	—	—	—	—	—	—	—	31.0	15.4	0.62	33.1
70	—	—	—	—	—	—	—	—	—	25.6	11.4	0.21	20.7
80	—	—	—	—	—	—	—	—	—	19.1	6.90	—	9.00
90	—	—	—	—	—	—	—	—	—	12.6	3.60	—	2.70
600	—	—	—	—	—	—	—	—	—	7.78	1.41	—	0.40
10	—	—	—	—	—	—	—	—	—	4.17	0.40	—	—
20	—	—	—	—	—	—	—	—	—	2.34	0.15	—	—
30	—	—	—	—	—	—	—	—	—	1.38	—	—	—
40	—	—	—	—	—	—	—	—	—	0.80	—	—	—
50	—	—	—	—	—	—	—	—	—	0.54	—	—	—
60	—	—	—	—	—	—	—	—	—	0.36	—	—	—
70	—	—	—	—	—	—	—	—	—	0.27	—	—	—
80	—	—	—	—	—	—	—	—	—	0.23	—	—	0.66
90	0.20	0.25	—	—	—	—	—	—	—	0.19	—	—	6.90
700	2.24	0.85	—	—	—	—	—	—	—	0.17	—	—	27.8
Luminous transmit.	2.8	2.4	2.8	0.78	1.86	0.88	0.69	0.36	0.26	20.1	9.0	0.032	31.4
Dominant wave lgth.	477.6	470.4	463.7	479.8	466.5	458.0	457.9	455.5	455.9	553.3	551.1	546.1	530.2
Excitation purity.	89.7	94.9	95.8	69.1	96.1	98.3	98.9	99.3	99.4	77.3	89.7	97.0	68.4

Wave length	Percent transmittance													
	No. 56	No. 57	No. 57A	No. 58	No. 59	No. 59A	No. 60	No. 61	No. 64	No. 65	No. 65A	No. 66	No. 67A	
400	—	—	—	—	—	—	—	—	9.00	—	—	—	12.3	1.10
10	—	—	—	—	—	—	—	—	9.20	—	—	—	13.0	0.93
20	—	—	—	—	—	—	—	—	8.75	0.23	—	—	15.0	1.28
30	—	—	—	—	—	0.16	—	—	9.20	0.61	0.16	—	18.4	3.16
40	—	—	0.19	—	—	0.37	—	—	11.3	1.58	1.32	—	23.2	6.40
50	—	—	0.87	—	0.40	1.26	0.19	—	15.5	4.10	5.50	—	31.2	10.5
60	0.16	0.44	2.56	—	1.90	4.57	1.38	—	23.3	9.00	13.0	—	42.2	17.7
70	3.12	3.10	7.80	0.23	7.70	13.2	5.38	—	34.4	16.8	24.9	—	55.5	28.5
80	13.0	13.1	21.6	1.38	21.0	30.0	15.0	0.33	46.8	24.9	36.6	—	68.4	41.4
90	34.5	31.9	41.7	4.90	41.5	50.8	32.0	4.00	56.6	31.3	45.1	—	77.6	52.1
500	59.0	50.5	58.8	17.7	59.0	66.0	48.4	16.6	62.1	33.7	45.8	—	82.7	57.9
10	73.0	60.6	67.9	38.8	67.7	73.0	57.2	32.3	62.9	32.4	39.7	—	84.6	58.8
20	79.0	63.3	70.1	52.2	69.8	75.1	59.2	40.0	59.1	27.5	29.7	—	84.0	55.4
30	79.9	61.0	67.6	53.6	67.2	73.2	55.5	39.6	51.6	20.7	17.8	—	82.6	47.5
40	77.5	55.0	61.8	47.6	61.5	68.5	47.5	34.5	41.3	13.7	7.90	—	79.1	36.6
50	72.6	47.1	53.5	38.4	54.0	62.0	36.8	26.3	28.0	6.50	2.40	—	73.7	25.0
60	66.1	37.3	43.3	27.8	45.0	54.4	25.2	17.3	16.2	1.66	0.32	—	67.1	14.2
70	58.0	26.5	31.6	17.4	35.0	44.5	14.4	9.70	7.95	0.40	—	—	58.8	5.50
80	46.1	16.6	19.4	9.0	24.0	33.0	6.3	4.40	3.10	—	—	—	47.2	1.40
90	33.8	8.69	9.70	3.50	14.0	22.0	1.82	1.66	0.80	—	—	—	34.5	0.28
600	24.0	3.70	4.50	1.50	7.95	14.6	0.48	0.38	—	—	—	—	24.4	—
10	18.7	1.60	2.00	0.41	4.90	10.5	0.10	—	—	—	—	—	18.5	—
20	13.2	0.49	0.87	—	2.70	6.92	—	—	—	—	—	—	13.7	—
30	7.22	—	0.22	—	1.00	3.16	—	—	—	—	—	—	7.70	—
40	3.02	—	—	—	0.17	1.07	—	—	—	—	—	—	3.00	—
50	1.48	—	—	—	—	0.50	—	—	—	—	—	—	1.46	—
60	1.91	—	—	—	—	0.91	—	—	—	—	—	—	1.91	—
70	7.95	—	—	—	—	3.00	—	—	—	—	—	—	6.17	—
80	23.0	—	0.16	—	4.00	10.0	—	—	—	—	—	—	19.9	—
90	44.1	—	1.15	—	12.0	20.0	2.10	—	0.10	—	0.20	—	42.6	—
700	64.8	—	3.17	0.53	22.6	30.0	8.70	—	4.50	—	2.18	—	63.1	0.40
Luminous transmit.	52.8	32.5	37.2	23.7	38.7	45.8	26.1	16.8	25.0	9.6	9.8	58.3	22.4	
Dominant wave lgth.	552.3	536.4	534.0	540.2	538.3	541.4	525.7	536.8	497.3	496.6	492.7	512.3	499.8	
Excitation purity.	78.2	69.2	62.1	88.1	66.0	59.3	62.2	85.4	55.0	67.8	77.4	21.5	55.8	

* Some transmission below 400 mμ. Consult the manufacturer.

TRANSMISSION OF WRATTEN FILTERS (Continued)

Percent transmittance

Wave length	No. 70	No. 72B	No. 73	No. 74	No. 75	No. 76	No. 77	No. 77A	No. 78	No. 78AA	No. 78A	No. 78B
400	—	—	—	—	—	0.22	—	—	37.2	43.0	56.0	64.1
10	—	—	—	—	—	0.18	—	—	41.7	46.0	58.6	66.5
20	—	—	—	—	—	0.29	—	—	44.2	48.7	61.0	68.4
30	—	—	—	—	—	1.38	—	—	44.6	49.8	61.8	69.5
40	—	—	—	—	—	3.50	—	—	44.2	49.7	61.8	70.0
50	—	—	—	—	—	3.50	—	—	41.7	48.0	61.0	69.4
60	—	—	—	—	1.97	1.92	—	—	38.0	44.9	58.7	67.5
70	—	—	—	—	10.0	0.51	—	—	33.8	40.3	55.0	65.4
80	—	—	—	—	17.4	—	—	—	27.5	35.6	51.0	62.9
90	—	—	—	—	18.0	—	—	—	23.5	30.9	47.1	59.8
500	—	—	—	—	13.0	—	—	—	19.5	26.5	43.5	57.0
10	—	—	—	0.96	7.35	—	0.30	0.10	15.8	23.4	40.0	54.2
20	—	—	—	7.95	3.20	—	9.10	5.35	13.8	20.3	36.9	51.4
30	—	—	—	14.6	0.83	—	13.5	1.90	11.8	17.8	34.4	49.3
40	—	—	—	12.9	0.14	—	46.1	35.0	10.5	16.6	32.7	48.1
50	—	—	—	7.60	—	—	78.0	71.8	9.56	14.9	31.2	46.7
60	—	—	2.24	3.06	—	—	75.8	63.1	8.53	13.2	29.4	45.0
70	—	—	5.97	0.83	—	—	8.00	—	7.77	12.1	28.0	43.6
80	—	—	4.56	0.12	—	—	1.00	—	7.41	11.6	27.5	43.1
90	—	1.26	2.00	—	—	—	0.32	—	6.93	11.1	27.0	42.9
600	—	5.89	0.56	—	—	—	16.2	1.60	6.45	10.40	26.0	41.8
10	—	5.25	0.10	—	—	—	52.1	32.1	5.50	9.20	24.1	40.0
20	—	2.88	—	—	—	—	83.0	78.0	4.80	7.70	21.8	37.6
30	—	1.26	—	—	—	—	84.9	79.5	3.94	6.50	19.7	35.5
40	—	0.48	—	—	—	—	88.1	86.5	3.46	5.60	18.6	34.2
50	0.63	0.14	—	—	—	—	89.8	89.2	3.24	5.50	18.4	33.6
60	10.5	—	—	—	—	—	89.8	89.0	3.16	5.60	18.5	34.0
70	35.0	—	—	—	—	—	85.5	79.5	3.39	5.80	18.7	34.1
80	55.2	—	—	—	—	—	76.1	62.5	3.45	6.10	19.0	34.5
90	70.0	—	—	—	—	0.13	75.0	62.4	3.51	6.10	19.3	34.8
700	79.0	—	—	—	0.14	1.24	86.5	83.0	3.90	6.50	20.2	36.0
Luminous transmit.	0.31	0.74	1.3	4.0	1.9	0.046	32.3	25.5	10.7	15.8	31.6	46.7
Dominant wave lgth.	675.6	604.9	574.9	538.6	487.7	449.2	579.9	581.5	471.1	473.4	475.7	477.2
Excitation purity.	100.0	100.0	100.0	96.7	90.4	99.7	99.0	99.1	63.0	54.5	33.7	20.7

Percent transmittance

Wave length	No. 78C	No. 79	No. 80A	No. 81	No. 81A	No. 81B	No. 81C	No. 81D	No. 81EF	No. 82	No. 82A	No. 82B
400	74.9	24.0	67.6	77.7	65.1	55.1	46.1	38.2	30.7	83.0	80.1	76.7
10	76.6	26.0	73.1	78.1	65.9	55.8	46.6	38.4	31.5	83.7	80.8	78.0
20	77.9	29.0	76.8	79.0	67.6	57.7	49.0	41.0	34.3	84.6	81.6	79.2
30	78.9	31.0	7.77	80.5	70.2	61.0	52.5	45.5	38.6	85.1	82.2	79.7
40	79.4	32.2	76.5	81.9	72.8	64.5	57.2	50.0	43.2	85.4	82.4	79.7
50	79.5	32.7	73.0	83.0	74.8	67.2	60.5	53.9	47.4	85.4	82.4	79.2
60	79.3	31.4	69.0	83.7	76.0	69.1	63.0	56.5	50.2	85.0	81.7	78.0
70	78.6	28.8	63.6	84.3	77.1	70.6	64.2	58.1	52.0	84.6	80.7	76.3
80	77.8	25.6	57.6	84.6	77.8	71.3	65.0	59.0	53.0	84.0	79.3	74.4
90	76.7	22.2	51.3	84.9	78.3	71.8	65.7	60.0	54.0	83.3	78.0	72.1
500	75.5	19.3	45.2	85.3	78.6	72.6	66.4	60.8	55.4	82.6	76.6	70.2
10	74.2	16.8	39.4	85.4	79.0	72.9	66.5	61.1	56.2	82.0	75.3	68.3
20	73.0	14.2	34.2	85.5	79.5	73.2	67.0	61.6	57.0	81.4	74.0	66.5
30	72.1	12.7	30.0	86.0	80.4	74.5	68.8	62.5	59.5	81.0	73.1	65.5
40	71.5	11.8	27.1	86.5	81.5	76.0	71.0	66.1	62.7	80.8	72.7	65.0
50	70.7	11.0	24.8	86.8	82.3	77.0	72.0	67.3	64.5	80.6	72.4	64.5
60	69.8	9.76	23.5	87.0	82.6	77.6	72.5	68.0	65.3	80.4	71.8	63.8
70	69.0	8.81	22.6	87.1	82.7	77.8	72.7	68.3	65.8	80.2	71.5	63.2
80	68.8	8.50	22.6	87.1	82.8	78.0	73.0	68.5	66.0	80.2	71.5	63.2
90	68.6	8.29	23.2	87.4	83.1	78.2	74.0	69.5	66.5	80.3	71.7	63.4
600	68.0	7.56	23.7	87.6	84.0	79.1	75.6	72.0	68.1	80.2	71.5	63.0
10	66.7	6.45	23.2	88.1	85.0	81.0	78.5	75.0	71.6	79.3	70.3	61.5
20	65.0	5.13	21.0	88.8	86.1	83.1	80.8	78.0	74.7	78.4	68.5	59.0
30	63.8	4.17	18.2	89.2	87.0	84.2	82.1	79.8	77.0	77.5	66.9	56.9
40	63.0	3.47	15.8	89.4	87.4	85.1	83.0	80.8	78.4	76.8	65.5	55.0
50	62.7	3.16	14.5	89.5	87.7	85.6	83.5	81.5	79.2	76.5	64.8	54.1
60	63.0	3.09	13.8	89.8	88.0	86.0	84.1	82.1	80.1	76.2	64.6	53.7
70	63.3	3.16	13.4	90.0	88.2	86.5	84.8	83.0	80.9	76.1	64.5	53.7
80	63.4	3.16	12.7	90.1	88.5	87.0	85.5	83.7	81.8	76.1	64.4	53.5
90	63.6	3.16	11.7	90.3	89.0	87.5	86.1	84.6	82.9	76.2	64.2	53.4
700	65.0	3.31	11.5	90.5	89.2	88.0	86.8	85.5	84.0	77.1	64.6	54.1
Luminous transmit.	70.4	11.3	28.4	86.8	82.0	76.9	72.0	67.4	64.0	80.7	72.5	64.6
Dominant wave lgth.	479.8	474.8	471.7	576.7	577.5	577.8	577.4	579.5	579.0	477.5	476.6	475.6
Excitation purity.	6.8	52.8	45.9	2.9	6.0	8.7	10.7	14.7	19.0	3.0	6.3	10.2

* Some transmission below 400 m μ . Consult the manufacturer.

TRANSMISSION OF WRATTEN FILTERS (Continued)

Wave length	Percent transmittance				Wave length	Percent transmittance			
	No. 87	No. 87C	No. 88A	No. 89B		No. 87	No. 87C	No. 88A	No. 89B
700	—	—	—	11.2	30	74.1	17.8	84.7	88.8
10	—	—	—	32.4	40	77.7	28.2	85.5	89.0
20	—	—	—	57.6	50	81.4	41.0	86.1	89.2
30	—	—	7.4	69.1	60	84.0	53.8	86.6	89.4
40	0.10	—	32.8	77.6	70	85.4	61.6	87.2	89.6
50	2.19	—	56.3	83.1	80	86.8	69.2	87.5	89.8
60	7.95	—	69.2	85.0	90	87.8	74.1	87.8	89.9
70	17.4	—	74.2	86.1	900	88.4	78.5	88.0	90.0
80	31.6	—	77.6	87.0	10	88.8	81.5	88.2	90.1
90	43.7	—	79.7	87.7	20	89.1	83.6	88.4	90.2
800	53.8	0.32	81.4	88.1	30	89.1	85.1	88.6	90.3
10	61.7	3.20	82.6	88.4	40	89.1	86.0	88.8	90.4
20	69.2	8.90	83.7	88.6	50	89.1	87.0	89.0	90.5

* Some transmission below 400 m μ . Consult the manufacturer.